

# Which table is right for you?

Surgical C-Arm Table 840



Surgical C-Arm Table 846



Pain Management C-Arm Table 870



Urology C-Arm Table 800



Brachytherapy C-Arm Table 810



3D Imaging C-Arm Table 820



- 058-840 Table, Surgical C-Arm, Contour, 115 VAC
- 058-840-10 Table, Surgical C-Arm, Rectangle, 115 VAC
- 058-845 Table, Surgical C-Arm, Contour, 230 VAC
- 058-845-10 Table, Surgical C-Arm, Rectangle, 230 VAC

- 058-846 Table, Surgical C-Arm, Contour, 115 VAC
- 058-846-10 Table, Surgical C-Arm, Rectangle, 115 VAC
- 058-847 Table, Surgical C-Arm, Contour, 230 VAC
- 058-847-10 Table, Surgical C-Arm, Rectangle, 230 VAC

- 058-870 Table, Surgical C-Arm, Contour, 115 VAC
- 058-870-10 Table, Surgical C-Arm, Rectangle, 115 VAC
- 058-875 Table, Surgical C-Arm, Contour, 230 VAC
- 058-875-10 Table, Surgical C-Arm, Rectangle, 230 VAC

- 058-800 Table, Urology C-Arm, 800, 115 VAC
- 058-805 Table, Urology C-Arm, 800, 230 VAC

- 058-810 Table, Brachytherapy, 810, 115 VAC
- 058-815 Table, Brachytherapy, 810, 230 VAC

- 058-820 Table, 3D Imaging C-Arm, 115 VAC
- 058-825 Table, 3D Imaging C-Arm, 230 VAC

## Shared Specifications:

**Dimensions:**  
Overall: 97" l x 26" w (246 x 66 cm) with OR accessory rails

**Tabletop:** 97" l x 24" w (246 x 61 cm) – choice of classic contoured or rectangular design

**Radiolucent Area:** 75" l x 24" w (191 x 61 cm)

**Tabletop Material:** Carbon fiber with integral head section

**Mattress:** Seamless, 2" thick (5 cm)

**Attenuation:** 1.2 mm Aluminum equivalence

**Patient Restraints:** Three body straps

**Finish:** Stainless Steel

**Patient Capacity:** 500 lb (227 kg)

**Warranty:** Three years parts and labor

**Certifications:** ANSI/AAMI ES60601-1:2005 +C1 (2009) +A1 (2012) +A2 (2010); CAN/CSA C22.2 No. 60601-1:14; CAN/CSA-C22.2 NO. 60601-2-46:12; IEC 60601-1-2-46:2016.

CDRH registered and FDA Device Listed



840

846

870

Motions:	840	846	870
<b>Height Adjustable:</b> 33.5" to 43.5" (84 to 109 cm)	—	Height Adjustable: 32" to 42" (81 to 106 cm)	—
<b>Tabletop X Motion (head-to-toe float):</b> 35" (88.9 cm)	—	Tabletop X Motion (head-to-toe float): 35" (89 cm)	—
<b>Tabletop Y Motion (side-to-side float):</b> 10" (25.4 cm)	—	Tabletop Y Motion (side-to-side float): 10" (25.4 cm)	—
<b>Trendelenburg:</b> 0° to 20°	—	—	Trendelenburg: 0° to 20°
<b>Reverse Trendelenburg:</b> 0° to 20°	—	—	Reverse Trendelenburg: 0° to 20°
<b>Isocentric Lateral Roll:</b> 0° to ±15°	—	—	Isocentric Lateral Roll: 0° to ±15°

Controls:	840	846	870
<b>Hand Control:</b> Activates height, isocentric lateral roll and Trendelenburg motions	—	Hand Control: Activates height	Hand Control: Activates height, isocentric lateral roll and Trendelenburg motions
<b>Foot Control:</b> Activates all motions (has release for float)	—	—	Foot Control: Activates all motions
<b>Mushroom-Top Control:</b> Releases X-Y Tabletop motions	—	Mushroom-Top Control: Releases X-Y Tabletop motions	—

Accessory Rails:	840	846	870
Standard OR accessory rails 49.5" x 1.12" x .375" (124 x 2.86 x .95 cm) mounted near foot end of table	—	Standard OR accessory rails 49.5" x 1.12" x .375" (124 x 2.86 x .95 cm) mounted near foot end of table	—
<b>Weight:</b> 641 lb (291 kg)	—	—	—
<b>Shipping Weight:</b> 865 lb (392 kg)	—	Shipping Weight: 800 lb (363 kg)	Shipping Weight: 749 lb (340 kg)
<b>Power:</b> 115 VAC or 230 VAC and battery. Battery automatically charges when table is plugged in. An optional battery with separate charging unit provides continuous, cordless operation.	—	Power: 115 VAC or 230 VAC	Power: 115 VAC or 230 VAC

Catheter Tray Extension (removable):	840	846	870
18" l x 24" w (45.7 x 61 cm) non-radiolucent	—	Optional	Optional

## Shared Specifications:

**Motions:**  
**Tabletop Y Motion (side-to-side):** 10" (25.4 cm)

**Trendelenburg:** 0° to 20°

**Reverse Trendelenburg:** 0° to 20°

**Lateral Roll:** 0° to 20°

**Controls:**  
**Hand Control:** Activates all motions

**Foot Control:** Activates all motions

**Accessory Rails:** Standard OR accessory rails, 1.125" x .375" (2.86 x .95 cm).

**Power:** 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger.

**Warranty:** Three years parts and labor

**Certifications:** ANSI/AAMI ES60601-1:2005 +C1 (2009) +A1 (2012) +A2 (2010); CAN/CSA C22.2 No. 60601-1:14; CAN/CSA-C22.2 NO. 60601-2-46:12; IEC 60601-1-2-46:2016.

CDRH registered and FDA Device Listed



800 and 810

820

Dimensions:	800 and 810	820
84.5" l x 25.5" w (214.6 x 64.8 cm) with OR accessory rails	—	Dimensions: 103.75" l x 26" w (264 x 66 cm) with OR accessory rails
<b>Tabletop:</b> 53" l x 21" w (134.6 x 53.3 cm); 84" l (213.4 cm) with extension	—	<b>Tabletop:</b> 103.75" l x 24" w (264 x 61 cm) includes head and leg extension
<b>Lightweight Extension (for patient transfer only, not for imaging):</b> 26.5" l x 21" w (67.3 x 53.3 cm)	—	—
—	—	<b>Radiolucent Area:</b> <b>3D Imaging Area:</b> 22.5" l x 24" w (57.1 x 61 cm)
<b>Radiolucent Extension:</b> 31" l x 25.5" w (78.7 x 64.8 cm) with OR accessory rails	—	<b>Radiolucent Extension (Fibersin™, non-3D):</b> 29" l x 22.5" w (73.6 x 57.1 cm)
<b>Radiolucent Area:</b> 54" l x 17.5" w (137.2 x 44.5 cm)	—	<b>Total Radiolucent Area (Combined):</b> 51.5" l x 24" w (131 x 61 cm)
—	—	<b>Head-End Extension (non viewable):</b> 18.2" l x 24" w (46 x 51 cm)
<b>Tabletop Material:</b> Three Fibersin™ sections – each with individual 2" mattress	—	<b>Tabletop Material:</b> Carbon fiber
<b>Mattress:</b> Individual, seamless, 2" thick (5 cm), one for each section	—	<b>Mattress:</b> Vinyl covered, seamless, 2" thick (5 cm)

Motions:	800 and 810	820
<b>Height Adjustable:</b> 34" to 43.5" (86.4 to 110.5 cm)	—	<b>Motions:</b> <b>Height Adjustable:</b> 34" h to 44" h (86 to 112 cm)
<b>Tabletop X Motion (head-to-toe):</b> 20" (51 cm)	—	<b>Tabletop X Motion (head-to-toe):</b> 20" (51 cm)
<b>Fowler Positioning:</b> Adjusts from 0° to 80° with touch handle lock/unlock	—	—

Optional:	800 and 810	820
<b>X-Ray Film Tray:</b> Removable film cassette tray for 14" x 17" film with cassette mask for 10" x 12" film. Film cassette tray can be removed from either side of the table.	—	—
<b>X-Ray Grid:</b> Removable X-Ray grid, 103 line, 8:1 ratio	—	—
—	—	<b>Head-End:</b> 30" l (75 cm)
—	—	<b>Foot-End:</b> 7.75" l (20 cm)
<b>Patient Restraints:</b> Two body straps, two ankle straps	—	—
<b>Finish:</b> Powder coat	—	—
<b>Patient Capacity:</b> In urologic position: 500 lb (227 kg) In use with radiolucent extension: 400 lb (181.4 kg) Weight tested to four times the patient load rating	—	<b>Patient Capacity:</b> 500 lb (227 kg); weight tested to four times the patient load rating.
—	—	<b>Weight:</b> 542 lb (246 kg)
<b>Shipping Weight 800:</b> 903 lb (409.6 kg)	—	<b>Shipping Weight:</b> 700 lb (318.1 kg)
<b>Shipping Weight 810:</b> 877 lb (398.6 kg)	—	—