



MIRION
TECHNOLOGIES

Dose Drawing

Syringe Shields

*Dose Drawing Syringe Shield
(glass non-replaceable)*



*Dose Drawing Syringe Shield
with replaceable glass*



FEATURES

- Lightweight, easy-to-use design
- Accommodates standard and safety syringes
- Constructed with a polymer material that acts as a shock absorber, the shields are resistant to breakage.

DESCRIPTION

Dose Drawing Syringe Shields reduce hand exposure when drawing doses from a shielded vial. Constructed of optically clear high density lead glass, the barrel of the syringe shield offers complete 360° visibility while reducing radiation exposure from Tc-99m by more than 99%.

The shield features quick and smooth syringe insertion with an O-ring seal and anti-roll cap.

ORDERING INFORMATION

Dose Drawing Syringe Shields:

Accommodates luer and non-luer lock syringes

- **007-661** Syringe Shield, 3 cc
- **007-663** Syringe Shield, 5 cc and 6 cc
- **007-665** Syringe Shield, 10 cc and 12 cc

Dose Drawing Syringe Shields with Replaceable Glass:

Accommodates luer and non-luer lock syringes

- **007-691** Syringe Shield, 3 cc
- **007-693** Syringe Shield, 5 cc and 6 cc
- **007-695** Syringe Shield, 10 cc and 12 cc

Replacement:

- **127-691** Syringe Shield Replacement Glass, 3 cc
- **127-693** Syringe Shield Replacement Glass, 5 cc and 6 cc
- **127-695** Syringe Shield Replacement Glass, 10 cc and 12 cc

