



MIRION
TECHNOLOGIES

Syringe Carrier

Single, Multiple, Small, Large and Shields



ORDERING INFORMATION

Shielded Syringe Carrier, .125" lead:

- **001-181** Syringe Carrier, Small
- **001-179** Syringe Carrier, Large

Shielded Syringe Carrier, .25" lead:

- **001-182** Syringe Carrier, Small
- **001-180** Syringe Carrier, Large

SPECIFICATIONS

001-179 Shielded Syringe Carrier, Large

Dimensions:

- I.D.: 8.25 x 3 x 2.9 in (21 x 7.6 x 7.4 cm) (l x w x h)
- O.D.: 9.5 x 4.4 x 3.5 in (24 x 11.2 x 8.9 cm) (l x w x h)

Lead Shielding:

- Sides, top and bottom: .125" thick (.32 cm)
- Ends: .25" thick (.64 cm)

Weight:

- 11.3 lb (5.1 kg)

001-181 Shielded Syringe Carrier, Small

Dimensions:

- I.D.: 8 x 1.9 x 1.97 in (20.3 x 4.8 x 5 cm) (l x w x h)
- O.D.: 9.25 x 3.4 x 2.6 in (23.5 x 8.6 x 6.6 cm) (l x w x h)

Lead Shielding:

- Sides, top and bottom: .125" thick (.32 cm)
- Ends: .25" thick (.64 cm)

Weight: 7

- .5 lb (3.4 kg)

DESCRIPTION

Shielded Syringe Carriers reduce exposure while storing or transporting radioactive material. The overlapping lid design with snap-latch closure prevents streaming. There are two sizes and thicknesses of lead to choose from. The ends of the carriers are double thick to reduce the exposure from the ends of syringes.

001-182 Shielded Syringe Carrier, Small

Dimensions:

- I.D.: 7.5 x 1.7 x 1.7 in (19 x 4.3 x 4.3 cm) (l x w x h)
- O.D.: 9.25 x 3.4 x 2.6 in (23.5 x 8.6 x 6.6 cm) (l x w x h)

Lead Shielding:

- Sides, top and bottom: .25" thick (.64 cm)
- Ends: .5" thick (1.3 cm)

Weight:

- 11 lb (4.9 kg)

001-180 Shielded Syringe Carrier, Large

Dimensions:

- I.D.: 7.7 x 2.6 x 2.7 in (19.6 x 6.6 x 6.9 cm) (l x w x h)
- O.D.: 9.5 x 4.4 x 3.5 in (24 x 11.2 x 8.9 cm) (l x w x h)

Lead Shielding:

- Sides, top and bottom: .25" thick (.64 cm)
- Ends: .5 in thick (1.3 cm)

Weight:

- 17 lb (7.7 kg)

